


PROCEDURE ROOM SET UP

UROLIFT®

CLINICAL & MEDICAL AFFAIRS

- ☐ Set up room for cystoscopy procedure
- ☐ UroLift® System ancillary scope kit
- ☐ Retrieval kit
- ☐ UroLift System Devices
- ☐ Privacy sheet and/or drapes
- ☐ Supplies adequately stocked
- ☐ Assigned staff responsibilities
- ☐ Review anesthesia protocol
- ☐ Fluid management (Irrigation and collection)
- ☐ Trash management (Procedure collection and disposal)

UROLIFT® SYSTEM TREATMENT ROOM



ANCILLARY SCOPE KIT

2.9 mm, 0°, 36 cm Cystoscope [UL-Scope]

20F Sheath [UL-Sheath]

Visual Obturator [UL-VO]

RIGID RETRIEVAL KIT

4 mm, 0°, 30 cm Cystoscope [UL-Scope4-0]

Rigid Scissor [UL-SCI-R]

Rigid Grasper [UL-GRASP-R]

SCOPE AND RETRIEVAL KIT TRAYS



Kit Sterilization Trays
[Procedure Tray: UL-PKTRAY]
[Retrieval Tray: UL-RRKTRAY]